

FLOW SAFE ® CPAP



MCHENRY WESTERN LAKE COUNTY EMS

CPAP Device: Why change?



- **OLD: Oxy-PEEP**

- Does not allow
- Meds via nebulizer
- PEEP measurement
- Shortages/recalls
- Packaged in lg box
- Headstrap may be uncomfortable & challenging to use
- 3 settings

- **New: Flow-Safe**

- Easier to repeat NTG
- Nebulizer can be used
- Actual PEEP measured
- Leaks detected
- Smaller package
- Comfortable headstrap
- 1 setting

Just a reminder...how it works



- CPAP works by “splinting” the lungs with a constant pressure of air
- This reduces work of breathing
- In CHF, forces excess fluid out of alveoli & interstitial space back into the vasculature
- Also decreases venous return to the heart thereby lessening its workload
- Bottom line: Takes patients close to needing intubation & rapidly reverses their condition avoiding complications & death



Flow-SAFE™ Disposable CPAP System



Rx

LARGE ADULT

Mercury Medical

11300 - 49th Street North
Clearwater, FL USA 33762-4807
800-237-6418 or 727-573-0088
Fax: 727-573-6040
www.mercurymed.com

EC REP CE

0086
Scanlan Group B.V.
Postbus 75664
Schiphol-Tripport 1118 ZS
The Netherlands

REF 10-57009

LOT 11249009

Additional Languages

(NL)	Spiegelde CPAP-systeem (Grote volwassenen)	(NO)	engangs-system for konstanttryk luftveisventilasjon (Stor voksenstørrelse)
(FR)	Système CPAP jetable (Grand adulte)	(PL)	System CPAP jednorazowego użytku (Osoba dorosła większego rozmiaru)
(DE)	Einweg-CPAP-System (Große Erwachsene)	(PT)	Sistema de CPAP descartável (Adulto de constituição forte)
(EL)	Συστημα CPAP μιας χρήσης (Μεγάλο άτομο, 18 ετών+)	(ES)	Sistema desechable de CPAP (Adulto grande)
(IT)	Sistema CPAP monouso (Adulti, misura grande)	(SV)	engångssystem för CPAP (Stor vuxen)

Made in Malaysia

#83-900-0329

1/2011

Rev. 2



CPAP



- **No Change**

- SOP
- Indications
- Contraindications

- **Change**

- Device
- Procedural steps

CPAP




- **Indications**

- Heart Failure - acute
- Asthma/COPD - severe
- Drowning - near
- Flail chest (w/o pneumothorax)


- **Contraindications**

- Airway
 - Unable to obtain adequate seal
 - Facial anomalies/injury
- Breathing
 - Inadequate resp rate/effort
 - Pneumothorax
 - Penetrating chest trauma
- Circulation
 - SBP <90 mmHg / DBP <60
- Disability - Consciousness
 - Decreased LOC
 - Unable to follow commands
- GI
 - Aspiration risk
 - Gastric distention
 - Vomiting
- Pregnant

THIS IS HOW CPAP FEELS TO THE PATIENT!!



FLOW SAFE[®] CPAP



OXYGEN FLOW RATE = PEEP

- 15 L = approx 3-4 cm H₂O PEEP
- 20 L = approx 6-7 cm H₂O PEEP
- 25 L = approx 8.5-10 cm H₂O PEEP



TITRATE PEEP TO PATIENT'S WORK OF BREATHING



- **OLD PROCEDURE**

- PEEP started at 10cm H₂O
- Not titrated to work of breathing

- **NEW PROCEDURE**

- **START AT 15 L WHICH EQUALS 3-4 cm H₂O PEEP**
- Titrate up to patient response (decreased work of breathing)
- Maximum PEEP 25 L or 8.5 - 10 cm H₂O PEEP


Flow-Safe® CPAP --- PEEP



- Pressure manometer
between CPAP valve & face mask
measures actual PEEP being delivered
- Leak in system
(e.g., lack of tight face-mask seal) will affect PEEP




What is the PEEP?



What is the PEEP?




Assess PEEP as patient exhales



EXPLAIN PROCEDURE TO PATIENT



- **MAY REQUIRE EXTENSIVE COACHING.....**
- **If they are anxious, a little Versed may help.**
- **See appropriate SOP**



OPEN PACKAGE



- **SECURELY CONNECT MASK TO VALVE/TUBING**



OPEN PACKAGE



- **Yellow label displays PEEP settings**



O2 Flow



- Attach CPAP O2 tubing to regulator/flow-meter
- Begin O2 flow @ 15 L
- If needed, slowly increase O2 to desired O2 sat/PEEP (do not exceed 30 LPM)



UNDO 1 OR 2 OF THE QUICK RELEASE CLIP(S)



Hold, or have pt. hold mask snugly to face...
Good face-mask seal is critical





Tighten head straps using Velcro tabs



Adjust forehead pad flat on forehead



Adjust forehead pad flat on forehead



USE QUICK RELEASE TO ADMINISTER NITRO PER SOP



SOP: If systolic is 90 or above, give nitro every 5 minutes unlimited

IF YOU RECEIVE PERMISSION TO USE CPAP WITH SEVERE ASTHMA/COPD & PATIENT IS WHEEZING



The patient will need a neb treatment



CONNECT NEBULIZER BETWEEN
CPAP VALVE & MASK

Capnography can be used under CPAP



CPAP with MWLCEMS Neb

Adaptors



- Some nebulizers may require an adaptor (same adapter used for in-line- nebs)
- You will need this for Good Shepherd nebulizers



Flow-Safe® CPAP Device



- **PEEP depends on**

1. O2 flow rate
2. Face mask seal (leak)

- **Oxygen concentration**

depends on

1. O2 flow rate
2. Respiratory rate
3. Tidal volume

How long tanks will last at different flow rates



Minutes of Oxygen by Cylinder Size
Based on full 2200 PSI Cylinders

Flow/L/Min.	D Cylinder	E Cylinder	M Cylinder
5	70	123	703
6	58	102	598
8	44	77	498
10	35	61	374
12	29	51	299
15	23	41	199
20	16	29	175
25	14	23	140

CPAP Device Transition



- Attrition (use your old one, replace with a new one)
- Devices already on vehicles will not immediately be replaced
- Transition will occur w/ replacement devices
- When an old device used, it will be replaced w/ new device (as hospital stock allows)
- New devices available after 12-1-11

Documentation



- Document the oxygen flow rate
- Document the PEEP reading on the manometer
- This can be done in your narrative
- Remember to use the manometer reading when patient exhales.

Thanks!



***THANKS TO THE NWC EMS SYSTEM AND
KATHY KNOP FOR SOME OF THE CONTENT
OF THIS PRESENTATION***