McHenry Western Lake County EMS | McHenry County College Skill Proficiency Record

< ECG Monitoring / 12-lead Acquisition >

Provider Name:

See Illustration
Appendix B

Date: _____

Performance Standard Verbalizes appropriate scope and breadth of ECG monitoring for skill level [] EMR - None [] EMT - Electrode placement and transmission [] Medic - Electrode placement, interpretation, treatment and transmission	Satisfactory Test 1 Test 2 Test 3			Required Coaching Test 1 Test 2 Test 3			Unsatisfactory		
		[]		[]	[]			[]	[]
Explains procedure to patient and obtains consent	. []	[]	[]	[]	[]	[]	[]	[]	[]
Preps skin by removing oil, dirt and other contaminants that would prevent proper electrode adhesion	[]	[]	[]	[]	[]	[]	[]	[]	[]
*Explains different limb lead placement locations based on different teachings and manufacturer recommendations. [] Manufacturer recommends distal placement on limbs [] Practical application calls for proximal placement on limbs [] Limb lead placement on trunk is not appropriate	[]	[]	[]		[]	[]	[]	[]	[]
*Properly places limb leads	[]	[]	[]	[]	[]	[]	[]	[]	[]
Verifies monitor connection of all leads and ensures proper good quality tracings.	[]	[]	[]	[]	[]	[]	[]	[]	[]
*Verbalizes causes of poor tracings and troubleshooting techniques. [] Old/dried out electrodes [] Dirty/diaphoretic/hairy skin [] Cable displacement / improper connection to monitor [] Road vibration / Patient tremors	[]	[]	[]	[]	[]	[]	[]	[]	[]
Demonstrates ability to print and transmit tracing.	[]	[]	[]	[]	[]	[]	[]	[]	[]
*Verbalizes need for 12-lead acquisition and transmission [] C/o chest pain or suspicion of ACS [] Post ROSC [] As indicated per protocols	[]	[]	[]	[]	[]	[]	[]	[]	[]

Continued on next page.

McHenry Western Lake County EMS | McHenry County College Skill Proficiency Record

< ECG Monitoring / 12-lead Acquisition (cont.) >

Performance Standard *Properly places all leads	Satisfactory Test 1 Test 2 Test 3	Required			
*Properly places all leads	And the same of th	Coaching Test 1 Test 2 Test 3	Unsatisfactory Test 1 Test 2 Test 3		
*Properly places all leads [] V1 – 4 th intercostal space, just right of the sternum [] V2 – 4 th intercostal space, just left of the sternum [] V3 – Centered between V2 and V4 over 5 th rib [] V4 – 5 th intercostal space, mid-clavicular [] V5 – 5 th intercostal space, centered between V4 and V6 [] V6 – 5 th intercostal space, mid-axillary		[][][]	[] [] []		
Verifies monitor connection of all leads and ensures proper good quality tracings.	[][][]	[] [] []	[][][]		
*Demonstrates ability to print and transmit tracing.	[][][]	[][][]	[][][]		
Performance Validation: [] Pass [] Retest	>	[] Pass	[] Fail		